

BIOMERICA
MATERIAL SAFETY DATA SHEET

PRODUCT: PSA Rapid Test (OTC)

CODE #: 1081-1

Eye Contact: No component in this product represents a danger to exposure in this manner.

Inhalation: No component in this product represents a danger to exposure in this manner.

Ingestion: Clear airway if the entire test device is ingested and seek medical attention if the individual is not conscious. If the desiccant packet or the test buffer has been ingested, give emetics to cause vomiting. A stomach siphon may be applied. Avoid milk and seek medical attention.

Section 5 -	FIRE AND EXPLOSION DATA
--------------------	--------------------------------

Suitable Extinguishing Media: Not applicable.

Extinguishing Media Not to be Used: None known.

Special Fire Fighting Procedures: None known.

Unusual Fire and Explosion Hazards: None known.

NFPA Rating: Not applicable.

Section 6 -	MEASURES IN CASE OF ACCIDENTAL SPILL
--------------------	---

Body Precautions: Not applicable.

Environmental Precautions: Not applicable.

Cleaning Measures: Not applicable.

Section 7 -	HANDLING AND STORAGE
--------------------	-----------------------------

Handling: No special precautions required.

Storage: Store kit reagents in 2 - 30° C.

Section 8 -	EXPOSURE CONTROL AND PERSONAL PROTECTION
--------------------	---

Body Protection: Not applicable.

Respiratory protection: Not applicable.

Hand Protection: Not applicable.

Eyes Protection: Not applicable.

Others: None known.

Section 9 -	PHYSICAL AND CHEMICAL CHARACTERISTICS
--------------------	--

**BIOMERICA
MATERIAL SAFETY DATA SHEET**

PRODUCT: PSA Rapid Test (OTC)

CODE #: 1081-1

Physical State: Solid Nitrocellulose Membrane contained in a plastic housing.

Color: White

Odor: None detectable.

pH: Not applicable.

Flash Point: Not applicable.

Flammability: Not Flammable.

Auto flammability: Will not occur.

Vapor Pressure: Not applicable.

Section 10 -	STABILITY AND REACTIVITY DATA
---------------------	--------------------------------------

Stability: Stable at room temperature.

Conditions to Avoid: None known.

Materials to Avoid: Not applicable.

Hazardous Decomposition Products: None Known.

Section 11 -	TOXICOLOGY INFORMATION
---------------------	-------------------------------

Effects of Dangerous Substances at the Concentrations Present in the Preparation:

Toxic Effects: None known.

Harmful Effects: None known.

Sensitizing Effects: None known.

Carcinogenic/Mutagenic/Teratogen Effects: None known.

Section 12 -	ECOLOGICAL INFORMATION
---------------------	-------------------------------

None Known.

Section 13 -	DISPOSAL CONSIDERATIONS
---------------------	--------------------------------

Disposal of Preparation Left: Dispose according to current regional and national rules.

Contaminated Containers: Dispose of following current regulations.

Section 14 -	TRANSPORTATION INFORMATION
---------------------	-----------------------------------

**BIOMERICA
MATERIAL SAFETY DATA SHEET**

PRODUCT: PSA Rapid Test (OTC)

CODE #: 1081-1

Road / Railway Haulage ADR/RID: Not restricted.
Sea Freight IMO (IMDG): Not restricted.
Air Freight IATA (ICAO): Not restricted.
UN Number: Not applicable.

Section 15 -	OTHER REGULATORY INFORMATION
---------------------	-------------------------------------

Classification: This preparation is classified and labeled according to EEC Directive n. 88/379.

Symbols of Danger: None.

"R" Phrases Indicating Specific Risks: None.

"S" Phrases Indicating Caution: None.

Section 16 -	OTHER INFORMATION
---------------------	--------------------------

Recommended Use: For research or *in vitro* diagnostic use only.

The information contained herein is believed to be correct as of the date of issue but is provided without warranty of any kind. Biomerica shall not be held liable for any damage resulting from handling or from contact with the above product. The recipient is responsible for observing any applicable laws and guidelines. Use the product according to the package insert.