



Safety Data Sheet

Issue Date: 06-Apr-2011

Revision Date: 16-Feb-2015

Version 3

1. IDENTIFICATION

Product Identifier

Product Name IgG Foods ELISA Kits

Other means of identification

SDS # BIO-011

Product Code

IgG 14 Foods Additives ELISA PRODUCT CODE #: 7151
IgG 90 Foods ELISA PRODUCT CODE#: 7180
IgG 45 Foods French ELISA PRODUCT CODE #: 7185
IgG 4 Foods ELISA PRODUCT CODE #: 7187
IgG 7 Foods ELISA PRODUCT CODE #: 7189
IgG 90 Foods ELISA PRODUCT CODE #: 7190
IgG 90 Med Foods ELISA PRODUCT CODE #: 7191
IgG 90 Custom Foods ELISA PRODUCT CODE #: 7193
IgG 14 Custom Foods ELISA PRODUCT CODE #: 7194
IgG 45 Custom Foods ELISA PRODUCT CODE #: 7195
IgG 21 Custom Foods ELISA PRODUCT CODE #: 7198
IgG 90 Foods Korean ELISA PRODUCT CODE #: 7199

Recommended use of the chemical and restrictions on use

Recommended Use For research or in vitro diagnostic use only.

Uses Advised Against Not for use on or in humans.

Details of the supplier of the safety data sheet

Supplier Address

Biomerica
17571 Von Karman Avenue
Irvine, CA 92614, USA

Emergency Telephone Number

Company Phone Number Phone: 1-800-854-3002
Email: bmra@biomerica.com
Emergency Telephone (24 hr) INFOTRAC 1-352-323-3500 (International)
1-800-535-5053 (North America)

2. HAZARDS IDENTIFICATION

Appearance Clear for Calibrators and
Green for Serum Diluent liquid
Red for Conjugate

Physical State Liquid

Odor None

Classification

Skin corrosion/irritation	Category 2
Serious eye damage/eye irritation	Category 1

Signal Word

Danger

Hazard Statements

Causes skin irritation
Causes serious eye damage



Precautionary Statements - Prevention

Use personal protective equipment as required
Wash face, hands and any exposed skin thoroughly after handling

Precautionary Statements - Response

If exposed or concerned: Get medical advice/attention
IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
Immediately call a poison center or doctor/physician
IF ON SKIN: Wash with plenty of soap and water
If skin irritation occurs: Get medical advice/attention
Take off contaminated clothing and wash it before reuse

Other Hazards

Calibrators are formulated with a buffer base, animal or human serum. The kit components that are made with human serum are tested by a United States Food and Drug Administration (USFDA) licensed method and found to be non-reactive for HIV-1, HIV-2, Hepatitis B surface antigen, and HCV. Because no test method can offer absolute assurance that these agents are absent, reagents should be handled at the Biosafety Level 2, as recommended for any potentially infectious human blood product, in the United States Center for Disease Control (USCDC) and National Institute of Health (USNIH) manual "Biosafety in Microbiological Laboratories", 1988. All bovine serum products used are derived from animals of US origin, processed in USDA licensed facilities
Wash Buffer, Serum Diluent and Calibrators in kits are formulated with Sodium Azide as a preservative. Concentrated Sodium Azide may react with copper and lead plumbing to form explosive metal azides. It may also react with acids to form explosive hydrazoic acid. If drain disposed, flush with large amounts of water to prevent azide build-up. Avoid contact with Stop Solution containing 1N H₂SO₄ (Sulfuric Acid), an irritant to the skin and mucous membranes. Substrate Solution A contains Dimethyl Sulfoxide. In case of contact with any of these reagents, wash thoroughly with water

3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	CAS No	Weight-%
Sulfuric Acid	7664-93-9	<5

If Chemical Name/CAS No is "proprietary" and/or Weight-% is listed as a range, the specific chemical identity and/or percentage of composition has been withheld as a trade secret.

4. FIRST-AID MEASURES

First Aid Measures

Eye Contact	Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Get immediate medical advice/attention.
Skin Contact	Take off contaminated clothing. Wash with soap and water.
Inhalation	Remove to fresh air. If person is having difficulty breathing, give oxygen and call a physician immediately.
Ingestion	Flush mouth with copious amounts of water, provided that the person is conscious, and seek medical attention.

Most important symptoms and effects

Symptoms Causes skin irritation. Causes serious eye damage.

Indication of any immediate medical attention and special treatment needed

Notes to Physician Treat symptomatically.

5. FIRE-FIGHTING MEASURES

Suitable Extinguishing Media

Chemical or water fire extinguisher.

Unsuitable Extinguishing Media None known.

Specific Hazards Arising from the Chemical

Calibrators, Serum Diluent and Wash Buffer in kits are formulated with Sodium Azide as a preservative. Sodium Azide may react with copper and lead plumbing to form explosive metal azides that are sensitive to mechanical shock, concussion, friction and sparks. May react with acids to form explosive hydrazoic acid. Dimethyl Sulfoxide in Substrate Solution A is combustible.

Hazardous Combustion Products None known.

Protective equipment and precautions for firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

Personal Precautions Wear gloves, impermeable shoe covers, and laboratory coat. Take care not to contaminate body. Ensure adequate ventilation.

Environmental Precautions Contain the spill to the smallest area possible. See Section 12 for additional Ecological Information.

Methods and material for containment and cleaning up

Methods for Containment Prevent further leakage or spillage if safe to do so.

Methods for Clean-Up Neutralize a Stop Solution spill with dilute base, then absorb the material with disposable towels. Soak a calibrator or control spill area with a 10% bleach solution and wipe up with disposable towels. Dispose of all contaminated trash in accordance with local regulations.

7. HANDLING AND STORAGE

Precautions for safe handling

Advice on Safe Handling Use personal protection recommended in Section 8. Take care not to splash, spill, or splatter standards, stop solution, or controls. Avoid contact with skin and eyes. Wash face, hands, and any exposed skin thoroughly after handling.

Conditions for safe storage, including any incompatibilities

Storage Conditions Store kit reagents in 2-8°C in refrigerators designated and labeled to contain human blood products.

Incompatible Materials Lead. Copper.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Exposure Guidelines

Chemical Name	ACGIH TLV	OSHA PEL	NIOSH IDLH
Sulfuric Acid 7664-93-9	TWA: 0.2 mg/m ³ thoracic fraction	TWA: 1 mg/m ³ (vacated) TWA: 1 mg/m ³	IDLH: 15 mg/m ³ TWA: 1 mg/m ³
Sodium Azide 26628-22-8	Ceiling: 0.29 mg/m ³ NaN ₃ Ceiling: 0.11 ppm Hydrazoic acid vapor	(vacated) S* (vacated) Ceiling: 0.1 ppm HN ₃ (vacated) Ceiling: 0.3 mg/m ³ NaN ₃	Ceiling: 0.1 ppm HN ₃ Ceiling: 0.3 mg/m ³ NaN ₃

Appropriate engineering controls

Engineering Controls Apply technical measures to comply with the occupational exposure limits.

Individual protection measures, such as personal protective equipment

Eye/Face Protection Wear approved safety goggles where a splash hazard exists.

Skin and Body Protection Wear suitable protective clothing.

Respiratory Protection Ensure adequate ventilation, especially in confined areas.

General Hygiene Considerations Handle in accordance with good industrial hygiene and safety practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Information on basic physical and chemical properties

Physical State	Liquid	Odor	None
Appearance	Clear for Calibrators and Green for Serum Diluent liquid		
	Red for Conjugate	Odor Threshold	Not determined
Color	Clear for Calibrators and Green for Serum Diluent		
	Red for Conjugate		

<u>Property</u>	<u>Values</u>	<u>Remarks • Method</u>
pH	5.5-7.4	Stop Solution pH <2
Melting Point/Freezing Point	Not applicable	
Boiling Point/Boiling Range	Not applicable	
Flash Point	Not applicable	
Evaporation Rate	Not determined	
Flammability (Solid, Gas)	Not flammable	
Upper Flammability Limits	Not applicable	
Lower Flammability Limit	Not applicable	
Vapor Pressure	Not determined	
Vapor Density	Not determined	
Specific Gravity	Not determined	
Water Solubility	100% soluble	
Solubility in other solvents	Not determined	
Partition Coefficient	Not determined	
Auto-ignition Temperature	Will not occur	
Decomposition Temperature	Not determined	
Kinematic Viscosity	Not determined	
Dynamic Viscosity	Not determined	
Explosive Properties	Not determined	
Oxidizing Properties	Not determined	

10. STABILITY AND REACTIVITY

Reactivity

Not reactive under normal conditions.

Chemical Stability

Stable under recommended storage conditions.

Possibility of Hazardous Reactions

None under normal processing.

Conditions to Avoid

None known.

Incompatible Materials

Lead. Copper.

Hazardous Decomposition Products

None known.

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure

Product Information	Potential biohazard
Eye Contact	Causes serious eye damage.
Skin Contact	Causes skin irritation.
Inhalation	Do not inhale.
Ingestion	Do not ingest.

Component Information

Chemical Name	Oral LD50	Dermal LD50	Inhalation LC50
Dimethyl sulfoxide 67-68-5	= 14500 mg/kg (Rat)	= 40 g/kg (Rat)	-
Sulfuric Acid 7664-93-9	= 2140 mg/kg (Rat)	-	= 510 mg/m ³ (Rat) 2 h
Polyoxyethylene Mono (Octylphenyl) Ether 9002-93-1	= 1800 mg/kg (Rat)	-	-
Sodium Azide 26628-22-8	= 27 mg/kg (Rat)	= 50 mg/kg (Rat) = 20 mg/kg (Rabbit)	-

Information on physical, chemical and toxicological effects

Symptoms Please see section 4 of this SDS for symptoms.

Delayed and immediate effects as well as chronic effects from short and long-term exposure

Skin corrosion/irritation Causes skin irritation.

Carcinogenicity IARC has classified "strong inorganic acid mist containing sulfuric acid" as a Category 1 carcinogen, substance that is carcinogenic to humans. This classification does not apply to liquid forms of sulfuric acid. Inorganic mist is not generated under normal use of this product.

Chemical Name	ACGIH	IARC	NTP	OSHA
Sulfuric Acid 7664-93-9	A2	Group 1	Known	X

Legend

- ACGIH (American Conference of Governmental Industrial Hygienists)
- A2 - Suspected Human Carcinogen
- IARC (International Agency for Research on Cancer)
- Group 1 - Carcinogenic to Humans
- NTP (National Toxicology Program)
- Known - Known Carcinogen
- OSHA (Occupational Safety and Health Administration of the US Department of Labor)
- X - Present

Numerical measures of toxicity

Not determined

12. ECOLOGICAL INFORMATION

Ecotoxicity

The product is not classified as environmentally hazardous. However, this does not exclude the possibility that large or frequent spills can have a harmful or damaging effect on the environment.

Component Information

Chemical Name	Algae/aquatic plants	Fish	Toxicity to microorganisms	Crustacea
Dimethyl sulfoxide 67-68-5	12350 - 25500: 96 h Skeletonema costatum mg/L EC50	34000: 96 h Pimephales promelas mg/L LC50 33 - 37: 96 h Oncorhynchus mykiss g/L LC50 static 40: 96 h Lepomis macrochirus g/L LC50 static 41.7: 96 h Cyprinus carpio g/L LC50		7000: 24 h Daphnia species mg/L EC50
Sulfuric Acid 7664-93-9		500: 96 h Brachydanio rerio mg/L LC50 static		29: 24 h Daphnia magna mg/L EC50
Sodium Azide 26628-22-8		0.8: 96 h Oncorhynchus mykiss mg/L LC50 0.7: 96 h Lepomis macrochirus mg/L LC50 5.46: 96 h Pimephales promelas mg/L LC50 flow- through		

Persistence/Degradability

Not determined.

Bioaccumulation

.

Mobility

Chemical Name	Partition Coefficient
Dimethyl sulfoxide 67-68-5	-2.03

Other Adverse Effects

Not determined

13. DISPOSAL CONSIDERATIONS

Waste Treatment Methods

Disposal of Wastes

Disposal should be in accordance with applicable regional, national and local laws and regulations.

Contaminated Packaging Disposal should be in accordance with applicable regional, national and local laws and regulations.

California Hazardous Waste Status

Chemical Name	California Hazardous Waste Status
Sulfuric Acid 7664-93-9	Toxic Corrosive
Sodium Azide 26628-22-8	Ignitable Reactive

14. TRANSPORT INFORMATION

Note Please see current shipping paper for most up to date shipping information, including exemptions and special circumstances.

DOT Not regulated

IATA Not regulated

IMDG Not regulated

15. REGULATORY INFORMATION

International Inventories

Chemical Name	TSCA	DSL	NDSL	EINECS	ELINCS	ENCS	IECSC	KECL	PICCS	AICS
Sulfuric Acid	Present	X		Present		Present	X	Present	X	X

Legend:

- TSCA - United States Toxic Substances Control Act Section 8(b) Inventory*
- DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List*
- EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances*
- ENCS - Japan Existing and New Chemical Substances*
- IECSC - China Inventory of Existing Chemical Substances*
- KECL - Korean Existing and Evaluated Chemical Substances*
- PICCS - Philippines Inventory of Chemicals and Chemical Substances*
- AICS - Australian Inventory of Chemical Substances*

US Federal Regulations

CERCLA

This material, as supplied, contains one or more substances regulated as a hazardous substance under the Comprehensive Environmental Response Compensation and Liability Act (CERCLA) (40 CFR 302)

Chemical Name	Hazardous Substances RQs	CERCLA/SARA RQ	Reportable Quantity (RQ)
Sulfuric Acid 7664-93-9	1000 lb	1000 lb	RQ 1000 lb final RQ RQ 454 kg final RQ
Sodium Azide 26628-22-8	1000 lb	1000 lb	RQ 1000 lb final RQ RQ 454 kg final RQ

SARA 313

Section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 (SARA). This product contains a chemical or chemicals which are subject to the reporting requirements of the Act and Title 40 of the Code of Federal Regulations, Part 372

Chemical Name	CAS No	Weight-%	SARA 313 - Threshold Values %
Sulfuric Acid - 7664-93-9	7664-93-9	<5	1.0
Sodium Azide - 26628-22-8	26628-22-8	<0.1	1.0

CWA (Clean Water Act)

This product contains the following substances which are regulated pollutants pursuant to the Clean Water Act (40 CFR 122.21 and 40 CFR 122.42)

Chemical Name	CWA - Reportable Quantities	CWA - Toxic Pollutants	CWA - Priority Pollutants	CWA - Hazardous Substances
Sulfuric Acid	1000 lb			X

US State Regulations

California Proposition 65

This product contains the following Proposition 65 chemicals.

Chemical Name	California Proposition 65
Sulfuric Acid - 7664-93-9	Carcinogen

U.S. State Right-to-Know Regulations

Chemical Name	New Jersey	Massachusetts	Pennsylvania
Dimethyl Sulfoxide 67-68-5	X		
Sulfuric Acid 7664-93-9	X	X	X
Urea Peroxide 124-43-6	X		
Sodium Azide 26628-22-8	X	X	X

16. OTHER INFORMATION

NFPA	Health Hazards Not determined	Flammability Not determined	Instability Not determined	Special Hazards Not determined
HMIS	Health Hazards Not determined	Flammability Not determined	Physical Hazards Not determined	Personal Protection Not determined

Issue Date: 06-Apr-2011
 Revision Date: 16-Feb-2015
 Revision Note: Updated Product List

Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

End of Safety Data Sheet